



Scientific / Metrology Instruments  
JEOL Mass Spectrometer Material analysis solution

Solutions for Innovation

# JEOL Mass Spectrometer

## Material analysis solutions



JEOL Ltd.

# JEOL offers solutions for your material analysis and applications

Our complete line of mass spectrometer systems and our expertise and years of experience help you solve the problems in analyzing functional polymer samples!



## Raw material examination

- Supplier selection
- Acceptance/



High sensitivity combined with high stability for simple, easy qualitative/quantitative analysis

### JMS-Q1500GC

Gas chromatography  
quadrupole mass spectrometer  
**GC-QMS**



Various ionization techniques and exact mass analysis for unknown material identification

### JMS-T200GC

Gas chromatography  
time-of-flight mass spectrometer  
**GC-TOFMS**



Simple yet robust all-round analysis

### JMS-T100LP

Liquid chromatography  
time-of-flight mass spectrometer  
**LC-TOFMS**



## Research and development



- Prototype testing
- Foreign object investigation
- Product examination

Material  
inspection

# JEOL Solutions


## Quality control



- Product complaints
- Problems with production process steps
- Final product evaluation

## Compliance

- WEEE/RoHS
- REACH
- European Directive 2005/84/EC



design for

mass spectrometer



Direct analysis of samples in various forms without preliminary treatment

### DART Ion Source

Ambient ionization  
time-of-flight mass spectrometer  
**DART-TOFMS**



High resolution observation of low to high molecular-weight ions

### JMS-S3000

Matrix assisted laser desorption/ionization  
time-of-flight mass spectrometer  
**MALDI-TOFMS**



## Prototype testing

GC-TOFMS estimates the composition and structure of known and unknown materials.


New materials used in prototypes often consist of compounds with unknown compositions and structures. Library searches by EI, which are designed for known materials, often fail to identify such compounds.

The JEOL GC-TOFMS, capable of 3 soft ionization techniques (FI, PI, and CI) and exact mass analysis, is ideal for the analysis of these prototype materials.

In this example, an unknown component in a liquid crystal mixture sample was analyzed.


A library search by EI failed to identify the component, as it was not present in the library. The GC-TOFMS, on the other hand, could identify the composition as well as the structure by utilizing molecular formula determination with CI and exact mass analysis of the EI fragment ions.

The GC-TOFMS makes it possible to determine the composition and structure of known and unknown materials.



EI and CI mass spectra of an unknown component in a liquid crystal mixture sample

Mode	Obs. $m/z$	Error (mDa)	Formula
EI <sup>+</sup>	69.0715	1.1	C <sub>5</sub> H <sub>9</sub>
	111.1187	1.3	C <sub>9</sub> H <sub>15</sub>
	195.0689	1.4	C <sub>13</sub> H <sub>9</sub> NO
CI <sup>+</sup>	333.1743	1.4	C <sub>22</sub> H <sub>23</sub> NO <sub>2</sub>
	334.1802	-0.5	C <sub>22</sub> H <sub>23</sub> NO <sub>2</sub>



Exact mass analysis results and estimated formula for an unknown component

## Product examination


GC-MS analysis combined with multiple classification analysis enables speedy identification of characteristic components.

While industrial polymer materials contain numerous compounds, those which “differentiate” the materials are limited. Multiple classification analysis, designed to survey multiple sets of data at a time, is effective in searching for such compounds.

In this example, 6 different vinyl acetate resins were analyzed by pyrolysis GC-EI-TOFMS.


PCA analysis of the data acquired separated the materials and extracted the characteristic compounds.

GC-MS analysis combined with multiple classification analysis enables speedy identification of characteristic components.



PCA score plot of 6 vinyl acetate resins

Characteristic compounds of sample B



PCA loading plot of 6 vinyl acetate resins

## Applications

- MSTips 040: The power of exact mass measurement: an example of unknown compound identification
- MSTips 108: Analysis of photo polymerization initiator in UV light curing adhesives by GC/TOFMS
- MSTips 220: Analysis of cyanoacrylate adhesive using the JMS-S3000 “SpiralTOF™”  
— Application of Kendrick Mass Defect plot analysis —



## Base material identification

MALDI-TOFMS enables precise structural analysis of synthetic polymers in polymer materials.


Physical properties of any polymer material change according to the structure of the oligomers, as well as the terminal groups, degree of polymerization, and molecular weight distribution. Analysis of the terminal group composition and degree of polymerization requires the direct observation of molecules.

MALDI-TOFMS, capable of observing low to high molecular weight samples with high resolution measurements, is ideal for this type of analysis.


In this example, a surfactant was analyzed.

The data shows 2 repeat structures, EO and PO, as well as their end structures and molecular weight distributions.

MALDI-TOFMS enables precise structural analysis of synthetic polymers in polymer materials.



MALDI mass spectrum of a surfactant



MALDI mass spectrum of a high molecular weight polymer

## Additives

Py-GC-MS analysis rapidly identifies additives in materials.


Today, a wide variety of additives are developed, produced, and distributed to provide polymers with new functions. Since these functions vary depending on the type and amount of additives in the polymer materials, it is critical to identify additives during material examination and quality control.

While pyrolysis (Py)-GC-MS is designed for direct analysis of materials which makes it effective in studying additives in polymers, this technique analyzes additives simultaneously with pyrolysates of matrix macromolecules. As a result, it can be difficult to identify all of these compounds by conventional one-dimensional chromatography. Comprehensive 2D GC (GCxGC) is effective for analyzing complex mixtures of components.

In this example, a commercial "X-ring" made of nitrile rubber (NBR) was analyzed.

The high resolution GCxGC measurement separated 4 additives from the numerous pyrolysate background components.

Py-GC-MS analysis combined with GCxGC makes it easier to identify additives.



Py-GCxGC-EI TIC chromatogram of nitrile rubber

## Applications

MALDI for Polymer Analysis: Synthetic Polymers and Additives

MSTips 138: Analysis of additives in plastic by thermal desorption (TD) GC-TOFMS: accurate mass measurement and isotope pattern matching

MSTips 111: Analysis of block copolymer by field desorption (FD) using JMS-T100GC "AccuTOF GC"



## Failure analysis

A variety of sampling techniques for GC-MS are effective in investigating product failures.

Quality control to investigate odor, foreign objects, and defects in products requires a number of analytical techniques since samples come in different types and forms.

GC-MS, a general purpose analytical system, supports multiple and diverse types of analysis and applications by combining different preliminary treatment devices.

In this example, gases produced from industrial rubber gloves were analyzed by HS-GC-MS using a head space (HS) technique for preliminary treatment. Four different gloves (two natural rubbers, a nitrile rubber, and a vinyl chloride resin) were heated at 50 °C, and the components of the resulting gases were compared. Each had a different profile with characteristic compounds detected from each glove sample.

A variety of preliminary sampling devices like HS are available for GC-MS, thus allowing the effective investigation of product failures.




Fig.3 TIC chromatograms of each samples@50°C


HS-GC-EI TIC chromatograms of 4 rubber gloves

## Surface deposit analysis

DART directly analyzes samples without preliminary treatment, and is effective in analyzing surface deposits and surface contamination.

DART (Direct Analysis in Real Time) is an ambient ionization technique capable of direct analysis of samples in different states such as gas, liquid, and solid.

In this example, cassette tapes were analyzed by DART to examine the surface deposits. Base peaks were detected at  $m/z$  199 for Tape A,  $m/z$  255 for Tape B, and  $m/z$  227 for Tape C. When these ions were subjected to exact mass analysis, the results suggested that these components were fatty acid used as a lubricant.



DART mass spectra of cassette tapes



Cassette tape being analyzed in DART

## Applications

Direct analysis of stains on cloth caused by cosmetic  
 MSTips 043: Analysis of Impurities in Ethylene glycols with JMS-T100GC AccuTOF GC



## RoHS



GC-MS is capable of speedy qualitative/quantitative analysis of substances on the RoHS list.

The polybrominated diphenyl ethers (PBDEs) and polybrominated biphenyls (PBBs) currently used as polymer flame retardants are suspected of being endocrine disruptors. The European Union specifically regulates the use of PBDEs and PBBs.


PBDEs and PBBs are bromine containing organic compounds that can be subjected to X-ray fluorescence (XRF) for primary screening to determine the concentration of total bromine. For resins and other materials found to have high concentrations during the primary screening, secondary screening with thermal extraction GC-MS and precision quantitative analysis with solvent extraction GC-MS are effective. Shown on the left is a flow chart for the analysis of PBDEs and PBBs.

In 2014, the RoHS directive was amended to include an additional 4 phthalate esters. Unlike PBDEs and PBBs, phthalate esters are organic compounds made up of C, H, and O so they cannot be subjected to primary screening with XRF. However, GC-MS is capable of analyzing phthalate esters, using the same techniques of simplified quantitative analysis and precise quantitative analysis that were used for PBDEs and PBBs.


The data shown here were acquired by thermal extraction GC-MS of a sample that contained the alternative materials DOTP, DIDP, and TOTM as plasticizers. Additionally, this sample was suspected to contain 4 additional phthalate ester impurities (DIBP, DBP, BBP, DEHP), which are on the amended RoHS list. This technique facilitates sampling with speedy qualitative/quantitative analysis at a rate of 30 minutes per sample.



Extracted ion chromatograms of a resin additive ( $m/z$  261: DOTP,  $m/z$  307: DIDP,  $m/z$  305: TOTM)



Chromatograms of 4 phthalate esters added to amended RoHS



Mass spectra and library search results of alternate materials (DOTP, DIDP, TOTM)



\*Appearance or specifications subject to change without notice.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.



3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577  
www.jeol.com ISO 9001 • ISO 14001 Certified

• **AUSTRALIA & NEW ZEALAND** /JEOL(AUSTRALASIA) Pty.Ltd, Suite 1, L2 18 Aquatic Drive - Frenchs Forest NSW 2086 Australia • **BELGIUM** /JEOL (EUROPE) B.V, Planet II, Gebouw B Leuvensesteenweg 542, B-1930 Zaventem Belgium • **BRAZIL** /JEOL Brasil Instrumentos Cientificos Ltda. Av. Jabaquara, 2958 5º andar conjunto 52 ; 04046-500 Sao Paulo, SP Brazil • **CANADA** /JEOL CANADA, INC. 3275 1ere Rue, Local #8 St-Hubert, QC J3Y-8Y6, Canada • **CHINA** /JEOL(BEIJING) CO., LTD. Zhongkeziyuan Building South Tower 2F, Zhongguancun Nansanjie Street No. 6, Haidian District, Beijing, P.R.China • **EGYPT** /JEOL SERVICE BUREAU 3rd Fl. Nile Center Bldg., Nawal Street, Dokki, (Cairo), Egypt • **FRANCE** /JEOL (EUROPE) SAS Espace Claude Monet, 1 Allée de Giverny 78290, Croissy-sur-Seine, France • **GERMANY** /JEOL (GERMANY) GmbH Gute Aenger 30 85356 Freising, Germany • **GREAT BRITAIN & IRELAND** /JEOL (U.K.) LTD, JEOL House, Silver Court, Watchmead, Welwyn Garden City, Herts AL7 1LT, U.K. • **ITALY** /JEOL (ITALIA) S.p.A. Palazzo Pacinotti - Milano 3 City, Via Ludovico il Moro, 6/A 20080 Basiglio(MI) Italy • **KOREA** /JEOL KOREA LTD, Dongwoo Bldg, 7F, 1443, Yangjae Daero, Gangdong-Gu, Seoul, 05355, Korea • **MALAYSIA** /JEOL(MALAYSIA) SDN.BHD, 508, Block A, Level 5, Kelana Business Center, 97, Jalan SS 7/2, Kelana Jaya, 47301 Petaling Jaya, Selangor, Malaysia • **MEXICO** /JEOL DE MEXICO S.A. DE C.V, Arkansas 11 Piso 2 Colonia Napoles Delegacion Benito Juarez, C.P. 03810 Mexico D.F., Mexico • **RUSSIA** /JEOL (RUS) LLC Krasnogorletarskaya Street, 16, Bld. 2, 127473, Moscow, Russian Federation • **SCANDINAVIA** /SWEDEN JEOL (Nordic) AB Hammarbacken 6A, Box 716, 191 27 Sollentuna Sweden • **SINGAPORE** /JEOL ASIA PTE.LTD, 2 Corporation Road #01-12 Corporation Place Singapore 618494 • **TAIWAN** /JIE DONG CO., LTD, 7F, 112, Chung Hsiao East Road, Section 1, Taipei, Taiwan 10023 (R.O.C.) • **THE NETHERLANDS** /JEOL (EUROPE) B.V. Lireweg 4, NL-2153 PH Nieuw-Vennep, The Netherlands • **USA** /JEOL USA, INC. 11 Dearborn Road, Peabody, MA 01960, U.S.A.

No.2101K760C(Bn)