

“Remainders of KM” plot for polymers using msRepeatFinder: compositional mapping over a broad mass range

Product used : MS

Mass Spectrometry (MS) with soft ionization such as matrix-assisted laser desorption ionization (MALDI) allows the compositional analysis of polymers (repeating units, chain terminations) of low dispersity. Combining a size exclusion chromatography (SEC) fractionation with a high-resolution MALDI SpiralTOF™ MS analysis enables the evaluation of the composition of polydisperse polymeric samples over a broad mass range (high-resolution/high-accuracy mass measurements in the low mass range <4 kDa, isotopic resolution in higher mass range < 30kDa). However, the higher the resolution the more the peaks detected in the mass spectrum with one spectrum per fraction, making the interpretation of mass spectral data the rate-limiting step of the whole analytical procedure. A “remainders of Kendrick mass” analysis (RKM) is proposed as a rapid post-acquisition data processing using visual maps from concatenated low/high-accuracy and low/high mass range data.

Experimental

A 1 mg mL⁻¹ solution of poly(ϵ -caprolactone) (PCL, Polymer Source, P1302-CL) in CHCl₃ was fractionated by SEC (HLC8220 GPC system, Tosoh, TSKgel multipore HXL-M columns, flow rate: 1 mL min⁻¹, 0.5 mL per fraction). Mass spectra were recorded with a JMS-S3000 SpiralTOF™ mass spectrometer (matrix: DCTB, no salt added). Plots were computed using msRepeatFinder 3.0.

SEC-MALDI-MS

Five fractions are collected from the SEC elution turning the PCL sample into aliquots of low dispersity (Fig. 1A). The main distributions in the last fraction (m/z 2000 to 5000, Fig. 1B) are assigned to sodiated (H, OH)-ended PCL (red circles) and (H, C₃H₇O)-ended PCL (blue triangle) taking full advantage of the high resolution of the analyzer. The SpiralTOF™ analyzer is still capable of producing mass spectra for the four other fractions of increasing molecular weight up to 20000 Da with a neat isotopic resolution (Fig. 2A). Mass spectra are concatenated using msRepeatFinder3.0 to display the full mass spectral data in one graph (Fig. 2B).




Fig. 1. (A) SEC chromatogram. (B) Mass spectrum of fraction 5.




Fig. 2. (A) Mass spectra of the five fractions. (B) Concatenated mass spectra.

Kendrick mass defect (KMD) plots

The regular KMD plot from the concatenated mass spectra barely separates the sodiated (H, OH)- and (H, C₃H₇O)-ended oligomers in the lowest mass range but the plot becomes unresolved while increasing the chain length (Fig. 3). In spite of the isotopic resolution reached by the SpiralTOF™ analyzer, the mass accuracy is not high enough for a regular KMD analysis. The resolution-enhanced KMD plot using a fractional base unit CL/113 (Fig. 4) successfully separates the four main distributions (sodiated and potassiated (H, OH)- and (H, C₃H₇O)-ended oligomers) over the whole mass range. The separating power is nevertheless decreasing while the molecular weight is increasing and the plot becomes fuzzy with a low quality of point alignment.





Fig. 3. Regular KMD plot [1] from the concatenated mass spectra of the five fractions (base unit: CL, C₆H₁₀O₂, 114.0681).





Fig. 4. Resolution-enhanced KMD plot [2] from the concatenated mass spectra of the five fractions using the “fraction base KMD” option (base unit: CL/113).



Reminders of KM (RKM) plot

The RKM plot of the concatenated mass spectra reveals a great compositional homogeneity of PCL throughout a broad 20 kDa mass range with (H, OH)- and (H, C₃H₇O)-ended chains (sodium and potassium adduction, Fig. 5A). An additional cyclic ion series is observed in the lowest mass range (fraction #5, blue square) typical of ring-opening and polycondensation synthesis routes. The detection of (H, ONa)-ended oligomers seen in the fraction #5 (violet square) further validates the presence of a COOH acidic end-group and confirms the (H, OH) assignment. As compared to the unresolved or fuzzy KMD plots, points are perfectly aligned in the RKM plot throughout the whole mass range in spite of using different external calibrants (one per fraction). The isotopic shift (¹²C → ¹³C) is clearly seen while increasing the chain length (from ¹²C to ¹³C₁₆ for the largest chains, Fig. 5B). The high separating power of the RKM plot allows a rapid filtering of a given ion series (full series or fixed isotopic composition) over the whole mass range using the “Grouping Mode” of msRepeatFinder (Fig. 6, red bars assigned to (H, OH)-ended PCL from 2 kDa to 20 kDa).




Fig. 5. RKM plots from the concatenated mass spectra. (A) Full plot with assignments of end-groups. (B) Detail of the sodiated (H, OH) ion series.





Fig. 6. Instant selection of the whole sodiated (H, OH)-PCL series throughout the 20kDa mass range (five fractions at once) using the “grouping mode” of msRepeatFinder.

Prospects

The RKM plots are compatible with high-accuracy and low-accuracy mass spectral data from the SpiralTOF™ analyzer (isotopic resolution) and the linear TOF analyzer (oligomeric resolution) of the JMS-S3000 mass spectrometer [3]. It is also compatible with multiple charging potentially observed in the MALDI-MS analysis of high molecular weight polymers [4].

Acknowledgment

This application note is written based on the results of a joint research project with Dr. Hiroaki Sato and Dr. Thierry Fouquet in Research Institute for Sustainable Chemistry, National Institute of Advanced Industrial Science and Technology (AIST).

References

- [1] H. Sato, S. Nakamura, K. Teramoto, T. Sato. *J. Am. Soc. Mass Spectrom.* **2014**, *25*, 1346–1355.
- [2] T. Fouquet, H. Sato. *Anal. Chem.* **2017**, *89*, 2682–2686.
- [3] T. Fouquet, T. Satoh, H. Sato. *Anal. Chem.* **2018**, *90*, 2404–2408.
- [4] T. Fouquet, R. B. Cody, Y. Ozeki, S. Kitagawa, H. Ohtani, H. Sato. *J. Am. Soc. Mass Spectrom.* **2018**, *29*, 1611–1626.

Copyright © 2018 JEOL Ltd.

Certain products in this brochure are controlled under the “Foreign Exchange and Foreign Trade Law” of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with “End-user’s Statement of Assurance” and “End-use Certificate” in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

JEOL **JEOL Ltd.**

3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577
www.jeol.com ISO 9001 • ISO 14001 Certified

• AUSTRALIA & NEW ZEALAND /JEOL(AUSTRALASIA) Pty.Ltd. Suite 1, L2 18 Aquatic Drive - Frenchs Forest NSW 2086 Australia • BELGIUM /JEOL (EUROPE) B.V. Planet II, Gebouw B Leuvensesteenweg 542, B-1830 Zaventem Belgium • BRAZIL /JEOL Brasil Instrumentos Cientificos Ltda. Av. Jabaquara, 2958 5º andar conjunto 52 ; 04046-500 Sao Paulo, SP Brazil • CANADA /JEOL CANADA, INC. 3275 1ere Rue, Local #8 St-Hubert, QC J3Y-8Y6, Canada • CHINA /JEOL(BEIJING) CO., LTD. Zhongkeziyuan Building South Tower 2F, Zhongguancun Nanshanjie Street No. 5, Haidian District, Beijing, P.R.China • EGYPT /JEOL SERVICE BUREAU 3rd Fl. Nile Center Bldg., Navai Street, Dokki, (Cairo), Egypt • FRANCE /JEOL (EUROPE) SAS Espace Claude Monet, 1 Allée de Giverny 78290, Croissy-sur-Seine, France • GERMANY /JEOL (GERMANY) GmbH Gute Aenger 3c 85356 Freising, Germany • GREAT BRITAIN & IRELAND /JEOL (U.K.) LTD. JEOL House, Silver Court, Watchmead, Welwyn Garden City, Herts AL7 1LT, U.K. • ITALY /JEOL (ITALIA) S.p.A. Palazzo Pacinotti - Milano 3 City, Via Ludovico il Moro, 6/A 20080 Basiglio(MI) Italy • KOREA /JEOL KOREA LTD. Dongwoo Bldg. 7F, 1443, Yangjae Daero, Gangdong-Gu, Seoul, 05355, Korea • MALAYSIA /JEOL(MALAYSIA) SDN.BHD. 508, Block A, Level 5, Kelana Business Center, 97, Jalan SS 7/2, Kelana Jaya, 47301 Petaling Jaya, Selangor, Malaysia • MEXICO /JEOL DE MEXICO S.A. DE C.V. Arkansas 11 Piso 2 Colonia Napoles Delegacion Benito Juarez, C.P. 03810 Mexico D.F., Mexico • RUSSIA /JEOL (RUS) LLC. Krasnoprolletarskaya Street, 16, Bld 2, 127473, Moscow, Russian Federation • SCANDINAVIA /SWEDEN JEOL (Nordic) AB Hammarbacken 6A, Box 716, 191 27 Sollentuna Sweden • SINGAPORE /JEOL ASIA PTE.LTD. 2 Corporation Road #01-12 Corporation Place Singapore 618494 • TAIWAN /JIE DONG CO., LTD. 7F, 112, Chung Hsiao East Road, Section 1, Taipei, Taiwan 10023 (R.O.C.) • THE NETHERLANDS /JEOL (EUROPE) B.V. Lireweg 4, NL-2153 PH Nieuw-Vennep, The Netherlands • USA /JEOL USA, INC. 11 Dearborn Road, Peabody, MA 01960, U.S.A.