

Delta Tips

NMDT_0029

Baseline Correction in 2D Spectra


NMR data processing software

Delta
NMR Software
v5.0



If baseline is distorted in a 2D spectrum, bands of noise may appear along F_1 - and/or F_2 -axis. Such noise complicates the analysis of low-intensity peaks and is a serious problem in 2D NOESY and ROESY spectra. The figures below demonstrate automatic baseline correction along F_2 -axis.

Select **Actions – Base Correct – Transpose 0** to apply automatic baseline correction in F_2 -domain. There are no parameters to set.



★ To apply baseline correction in F_1 -domain, select **Transpose 1**.