

Delta Tips

NMDT_0042

Data Slate (Part 5)

NMR data processing software

Delta
NMR Software
v5.0



< How to display multiple 2D spectra in one Data Slate window >

It is possible to overlay two or more 2D spectra in one Data Slate window. It is a very useful function if you need to closely compare several 2D spectra.

In the example bellow, we will demonstrate how to overlay HMQC (A) and HMBC (B) spectra. At first, open the Spectrum A in one Data Slate window and the Spectrum B in another one. Next, follow the steps as follows:



- ① Push the **Finger** button  and ② the **Overlay** button  in the Data Slate window of the Spectrum A. Note that the cursor has changed into the **Finger** symbol .
- ③ Click into the Spectrum B with the cursor .

The Spectrum B has been loaded and displayed over the Spectrum A.

Spectrum A



Spectrum B



- ★ If you wish to adjust contour levels of Spectrum A or Spectrum B, click the tab of the spectrum, and then adjust contours as usual.

 1: cahe_HMQC...jdf  2: cahe_HMBC...jdf

Refer to User's Manual - Data Processing User's Manual for detail.

- ★ It is possible to remove the spectrum by selecting it and pushing the Remove button .

- ★ For 1D spectra, refer to NMDT_0025 and NMDT_0026.