

Operando-analysis of organic photovoltaic devices via ESR and EDMR methods

Related product : Electron Spin Resonance (ESR)

Characteristics and issues in organic photovoltaic devices

Organic photovoltaic devices (OPVs, see Fig. 1(a)) are light, thin, and cost-effective. Furthermore, OPVs have many unique properties which conventional inorganic semiconductor cannot achieve. Therefore, they are attracting considerable attention as next-generation renewable energy devices^[1]. The most important issues that inhibit commercialization are durability and stability (see Fig. 1(b)). To explore the origin of degradation in OPVs, precise analysis is needed on the view points in molecular scope.

ESR and EDMR




Fig. 2 (a) Targets and principle of ESR measurement and (b) of EDMR

Operando-ESR and EDMR

The observed OPV using blended P3HT:PC₆₁BM film as an active layer improves the efficiency of electron transport due to the insertion of a thin LiF layer between the active layer and the Al electrode. It is well known that post-annealing of this device, which is an essential manufacturing process, induces significant performance degradation^[2]. As shown in Fig. 3(a), this ESR spectrum reveals PC₆₁BM anion radicals generated at the interface between PC₆₁BM and a LiF / Al electrode via a chemical reaction. It is considered that these anions act as charge-scattering centers. In contrast, the EDMR spectrum of the same device, as shown in Fig. 3(b), reveals a newly found current path through different paramagnetic species, which are identified as P3HT polarons acting as recombination centers. These two spectra imply that different spins contribute to different mechanisms of device performance degradation^[3, 4].




Fig. 3 (a) ESR spectrum of a post-annealed OPV (b) EDMR spectrum of the same OPV
Both spectra were measured under DC bias conditions.

References

- [1] L. Ding, Organic Solar Cells: Materials Design, Technology and Commercialization (Wiley, 2022). [2] F. Li, J. Zhao, K. Yao, and Y. Chen, Chem. Phys. Lett. **553**, 36 (2012). [3] D. Liu et al., Appl. Phys. Lett. **104**, 243903 (2014). [4] T. Suzuki and K. Marumoto, J. Appl. Phys. **135**, 075002 (2024).

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.