

Rapid analysis of phthalates for RoHs2 by Py-GC-MS

Related products: Mass spectrometer (MS)

Overview

Phthalates are used as plasticizers for resins, but their use is restricted by various regulations due to concerns about their effects on the human body as endocrine disruptors. The revised RoHS Directive (RoHS2), which came into effect on July 22, 2019, regulates four phthalats (diisobutyl phthalate: DIBP, dibutyl phthalate: DBP, butyl benzyl phthalate: BBP, di-2-ethylhexyl (DIBP, DBP, butyl benzyl phthalate: BBP, di-2-ethylhexyl phthalate: DEHP) were regulated. As measurement methods for these substances, the solvent extraction-GC/MS method for precise analysis and the Py-GC/MS method for screening are specified in IEC-62321-8. We have attempted a rapid analysis of Py-GC/MS method using JMS-Q1500GC by optimizing Py and GC oven conditions, and achieved a cycle time of 14.75 min per measurement. Here, we report the results.



Figure1. Cycle time details

Experiment

PVC (NMIJ CRM 8152-a, about 0.1%), a certified reference material for the analysis of phthalates, was used as the measurement sample. A 0.50 mg sample was measured by Py-GC/MS to confirm the peak separation of the 4 phthalates and the repeatability of the area values at n=5. Table 1 shows the measurement conditions for Py-GC/MS.

PY			MS : JMS-Q1500GC			
Sample weight	0.50mg		lon source temp.	250°C		
Furnace Temp.	$150^{\circ}C \rightarrow 200^{\circ}C/min \rightarrow 320^{\circ}C(0.5min)$		Interface temp.	300°C		
GC			Ionization	El, 70eV		
Column	ZB-1HT 15m length, 0.25mmi.d. 0.1µm film thickness (Phenomenex Inc.)		Measurement mode	Scan m/z 50 to 600 with SIM DIBP: m/z 149, 205, 223 DBP: m/z 149, 205, 223 BBP: m/z 149, 205, 223		
Gas flow	He 1.0mL/min, constant flow			DEHP: m/z 149 167 279		
Injection mode	Split 1/50			DEIN: 1102 110, 101, <u>Ero</u>		
Inlet Temp.	300°C			*Under line : Quantitative i		
Oven temp.	80°C→50°C/min→200°C					
	→25°C/min→350°C					

Table 1. PY-GC/MS measurement conditions

Results

The TIC chromatogram is shown in Figure 2. The 4 phthalates could be well separated.



Next, the extracted ion chromatogram (EIC) of the quantitation ion for each of the 4 phthalates is shown in Figure 3. All components showed good detection sensitivity.



Figure3. Extracted Ion Chromatograms of phthalates in PVC

Finally, the area values of each measured sample at n=5 are shown in Table 2. For each phthalate, good results were obtained with a coefficient of variation of less than 10%.

	#1	#2	#3	#4	#5	Average	Standar d deviatio n	CV%
DIBP	178051	168249	158209	157911	167761	166036	8359	5.0%
DBP	166170	150082	152769	142665	148792	152096	8695	5.7%
BBP	210573	212346	205639	205776	243147	215496	15735	7.3%
DEHP	110307	108873	97884	93277	102132	102495	7212	7.0%

Table 2. Repeatability of area value

Summary

JEO

JEOL Ltd.

A rapid analysis of four phthalates, which are regulated by RoHS2, by Py-GC/MS method was attempted and the data was obtained with a cycle time of 14.75 min based on good detection sensitivity and repeatability.

Copyright © 2022 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577 www.jeol.com ISO 9001 • ISO 14001 Certified



