

Figure 2 shows TICC enlargements and volcano plots for retention times from 5 to 13 minutes obtained by the difference analysis function of msFineAnalysis iQ. A number of minute peaks were detected overlapping with the broad peak of glycerol at a retention time of around 9 minutes. Table 2 shows the integrated analysis result list. From these results, 11 kinds of compounds were qualitatively obtained through both samples.

A volcano plot is a scatter plot that can visualize characteristic compounds between samples, with the intensity ratio (Log₂(B/A)) on the horizontal axis and the statistical reproducibility (-Log₁₀(p-value)) on the vertical axis respectively. This time, the left area of the volcano plot shows the compounds specifically included in cyan and the right area in magenta. The compounds specifically contained in each sample were caprolactam for magenta, and four compounds for cyan, including ID: 005 "Ethanol, 2,2'-oxybis- (diethylene glycol)".

Figure 3 shows the mass spectrum of the peak of [ID:007]. An ion of m/z 121, presumed to be a molecular ion, could be detected on both EI and PI mass spectra. Table 3 shows the integrated analysis result list (top 5 candidates) by msFineAnalysis iQ. From this result, the peak [ID:007] was presumed to be "Benzenamine, N-ethyl- ". This compound was presumed to be a kind of solvent contained in water-based ink.




Figure 2 Volcano plot of variance compound analysis result

Table 2 Integrated qualitative analysis result of peak

Chromatogram information					Difference analysis result			Integrated analysis result								
ID	RT [min]	RI [iu]	Area	Height	Height [%]	Type	Log ₂ (B/A)	p-value	Compound	CAS#	Similarity	Lib. RI [iu]	ΔRI [iu]	Formula	Mw	Isotope Matching
001	1.23	645	57.66	22234643	77.67	A = B	0.35	0.202	Water	7732-18-5	904	317	329	H ₂ O	18	0.99
002	1.30	650	22.43	11165988	39.01	A = B	-0.54	0.265	Ammonia	7664-41-7	749	N/A	N/A	H ₃ N	17	N/A
003	1.36	653	0.22	1070288	3.74	B Only	> 4	0.184	Ethanol	64-17-5	920	427	226	C ₂ H ₆ O	46	0.83
004	1.43	658	16.47	28626418	100.00	A = B	0.10	0.463	Isopropyl Alcohol	67-63-0	953	489	169	C ₃ H ₈ O	60	0.99
005	6.82	984	0.82	153743	0.54	A Only	< -4	0.020	Ethanol, 2,2'-oxybis-	111-46-6	934	927	57	C ₄ H ₁₀ O ₃	106	0.33
006	6.95	992	0.24	52891	0.18	A Only	< -4	0.246	Aniline	62-53-3	832	977	17	C ₆ H ₇ N	93	0.85
007	9.33	1135	0.22	197769	0.69	A Only	< -4	0.012	Benzenamine, N-ethyl-	103-69-5	962	1128	8	C ₈ H ₁₁ N	121	0.97
008	10.16	1189	26.42	24601179	85.94	A = B	-0.26	0.073	Ethanol, 1-(2-butoxyethoxy)-	54446-78-5	953	1187	2	C ₈ H ₁₈ O ₃	162	-
009	11.34	1270	0.67	236129	0.82	B Only	> 4	0.184	Caprolactam	105-60-2	955	1253	17	C ₆ H ₁₁ N O	113	0.93
010	13.23	1408	0.50	273662	0.96	A = B	0.13	0.431	2,4,7,9-Tetramethyl-5-decyn-4,7-diol	126-86-3	949	1407	1	C ₁₄ H ₂₆ O ₂	226	-
011	14.01	1469	0.28	106402	0.37	A Only	< -4	0.072	Ethanol, 2-[2-(2-butoxyethoxy)ethoxy]-	143-22-6	938	1486	17	C ₁₀ H ₂₂ O ₄	206	-




Figure 3 Mass spectra of peak [007]

Table 3 Integrated qualitative analysis result of peak [007]

#	Compound	CAS#	Similarity	Similarity Reverse	Lib. RI [iu]	ΔRI [iu]	Formula	EI Base Peak (Lib.)	Mw	Mw confirmation	Isotope Matching
★ L01	Benzenamine, N-ethyl-	103-69-5	962	962	1128	8	C ₈ H ₁₁ N	106	121	?	0.97
L02	Benzenamine, 2-ethyl-	578-54-1	870	870	1122	14	C ₈ H ₁₁ N	106	121	?	0.97
L07	m-Ethylaniline	587-02-0	803	803	1141	5	C ₈ H ₁₁ N	106	121	?	0.97
L09	Pyridine, 3-ethyl-5-methyl-	3999-78-8	794	795	645-1355	0	C ₈ H ₁₁ N	106	121	?	0.97
L11	2,6-Xylidine	87-62-7	749	749	1168	32	C ₈ H ₁₁ N	121	121	?	0.97

Conclusion

In this report, we reported an example of difference analysis by msFineAnalysis iQ for major constituents in two colors (cyan and magenta) of water-based inks. msFineAnalysis iQ uses not only library DB search but also multiple identification functions such as retention index and isotope matching, so highly accurate qualitative analysis is possible. This software is expected to improve qualitative accuracy and efficient analysis work in GC-QMS analysis.

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

