

## Analysis of volatile organic compounds inside a vehicle using thermal desorption-GC-MS

Product used: Mass Spectrometer (MS)

### Introduction

Thermal desorption (TD) is a GC pretreatment device for introducing gas samples such as atmospheric air and indoor air. It can also be used for gases generated from solid or liquid samples. Figure 1 shows an overview of TD operation. ① Sampling: Using a pump, etc., the gas sample is passed through and collected in a sample tube containing an adsorbent. ② Tube desorption: The sample tube is heated in the TD device, and the desorbed gas is collected in a cooled focusing trap. ③ Trap desorption: The focusing trap is heated at high speed and the desorbed gas is introduced into the GC in a narrow band width. These two steps of thermal desorption make it possible to obtain high sensitivity and sharp chromatogram peaks.

In this MSTips, we introduce the analysis of volatile organic compounds (VOCs) inside a vehicle, which is one of the main applications of TD-GC-MS. TD-100xr (Markes International Ltd) and JMS-Q1600GC (JEOL) were used in the experiment. The analysis conditions were based on the international standards ISO 12219-1 and JASO Z125.




Figure 1 Overview of TD operation

### Features of TD-100xr

- DiffLok cap: No need to remove from sample tube during measurement
- Electronically cooled focusing trap: Cools down to -30 °C without refrigerant
- Sample recollection: Collect the split flow during trap desorption and reuse the sample




Figure 2 Features of TD-100xr

## Target component


We measured 7 components of VOCs which have guideline values of concentration inside a vehicle and are specified TD-GC-MS as the measurement method. We also measured three components which are under additional consideration. Table 1 shows the guideline values of concentration established by the Ministry of Health, Labor and Welfare.

**Table 1 Guideline values of concentration ( Left : current 7 components, Right : additional 3 components)**

Name	Concentration ( $\mu\text{g}/\text{m}^3$ )	Name	Concentration ( $\mu\text{g}/\text{m}^3$ )
Toluene	260	2-Ethyl-1-hexanol	130
Xylene	200	Texanol	240
Ethylbenzene	3800	2,2,4-Trimethyl-1,3-pentanediol diisobutyrate : TXIB	100
Styrene	220		
Tetradecane	330		
Dibutyl phthalate : DBP	17		
Bis(2-ethylhexyl)phthalate) : DEHP	100		

## Experiment

For the current 7 components, 1, 10, and 100 ng of mixed standard solutions were added to sample tubes. For the additional 3 components, 50 ng of the mixed standard solution was added to the sample tube. The sampling volume of 3L described in ISO12219-1 corresponds to concentrations of 0.3, 3, 30  $\mu\text{g}/\text{m}^3$  and 15  $\mu\text{g}/\text{m}^3$ . Additionally, 47 ng of Toluene- $d_8$  was added to the sample tube as an internal standard. By using a special tool when adding these standard solutions, it is possible to easily collect the entire amount (Figure 3). Table 2 shows the measurement conditions of TD-GC-MS. In this measurement, we measured the current 7 components and additional 3 components separately, but these can also be measured simultaneously.



**Figure 3 Solution loading tool**

**Table 2 Measurement conditions**

TD conditions		GC conditions	
Thermal Desorption	TD-100xr (Markes International Ltd)	Gas Chromatograph	8890A GC (Agilent Technologies)
Sample tube type	Tenax TA	Column	ZB-5MSi (Phenomenex) 30m x 0.25mm, 0.25 $\mu\text{m}$
Tube desorption	280°C (10min), 30mL/min, Splitless	Oven Temperature	40°C(2min)-15°C/min -320°C(10min)
Focusing trap type	General purpose Hydrophobic (T2)	Carrier flow	He, 1.0mL/min
Trap cooling	0 °C	MS conditions	
Trap desorption	280°C (3min), 20mL/min, Split 20:1	Spectrometer	JMS-Q1600GC (JEOL Ltd.)
Flow path temperature	200 °C	Ionization	EI+:70eV, 50 $\mu\text{A}$
		Ion source temperature	250 °C
		Mass range	Scan, $m/z$ 45-500

## Measurement results of current 7 components

Figure 4 shows the TIC chromatogram at 100 ng/tube. Since xylene has 3 isomers as 2 peaks, the total number of peaks is 8.




Figure 4 TIC chromatogram of 100 ng/tube

Figure 5 shows the extracted ion chromatograms (EICs) and calibration curves of DBP ( $m/z$  149) and DEHP ( $m/z$  149). Good peak shapes and S/N sensitivity were obtained for all components at the lowest concentration. The determination coefficients of the calibration curve were 0.999 or higher, indicating good linearity.






Figure 5 EICs chromatograms and calibration curves (Left : DBP, Right : DEHP)

## Measurement results of additional 3 components

Figure 6 shows the TIC chromatogram of 50 ng/tube. Texanol has two isomers as two peaks, the total number of peaks is four.




Figure 6 TIC chromatogram of 50 ng/tube

Figure 7 shows the extracted ion chromatograms (EICs) and mass spectra of each component. Good peak shapes and S/N sensitivity were obtained for all components. \*Actually, Texanol-1 and Texanol-2 have different structures, but only one structure was registered in the library database, and both peaks matched that structure.




Figure 7 EICs and mass spectra

## Conclusion

Using the TD-100xr and JMS-Q1600GC, it is possible to analyze the current 7 components and additional 3 components of VOCs inside a vehicle with high sensitivity and linearity. These devices are expected to be useful in analyzing gas samples such as atmospheric air and indoor air.

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

