

Utilizing Python® in JASON via HDF5 Files and the External Command Function

Product used : Nuclear Magnetic Resonance (NMR)

Currently, JASON[1] provides an efficient environment for NMR spectral analysis using Python[2] through its hierarchical HDF5 file structure and integration with BeautifulJASON.

This document introduces the following two key topics:

- Overview of HDF5 files and their application in JASON
- Use cases of the External Command function

1. Overview of HDF5 Files and Their Application in JASON

What is HDF5?

HDF5 (Hierarchical Data Format version 5) is a file format designed for the efficient storage and management of large and complex datasets. Widely used in scientific and engineering fields, HDF5 offers the following key features:

- **Hierarchical Structure:** Data is organized into groups and datasets in a directory-like format, making it easy to manage large volumes of information.
- **Support for Large Datasets:** Optimized for storing and accessing massive data, HDF5 is ideal for scientific computing and image analysis.
- **Wide Language Support:** HDF5 can be accessed from a wide range of scientific programming languages, including C, C++, Fortran, Python, MATLAB, R, Java, Julia, and many others, thanks to extensive native and community-supported libraries.
- **Cross-Platform Compatibility:** The HDF5 format and libraries are fully cross-platform, with mature support for Windows, macOS, and Linux. HDF5 can also be built for mobile platforms such as iOS and Android for specialized applications, and browser-based access is possible through JavaScript/WebAssembly libraries.
- **Random Access:** Enables efficient read/write operations on specific data segments, making it suitable for large-scale data processing.

Advantages of HDF5 in NMR Data Analysis

- **Efficient NMR Spectral Analysis:** Enables fast and selective extraction of multidimensional data, facilitating smooth analysis workflows.
- **Hierarchical Data Access:** For example, data can be retrieved from paths such as `JasonDocument/NMR/NMRData/0/DataPoints` or `SpecInfo`.
- **Metadata Utilization:** Attributes stored in `SpecInfo[3]` (e.g., `SW`, `SpectrometerFrequencies`, `SpectrumRef`) can be easily accessed to obtain essential experimental conditions.

Table 1. Example Workflow for NMR Spectral Analysis

Step	Process Description	Libraries Used	Representative Code Snippet (Excerpt)
1. Data Loading	Retrieve NMR data from .j5 file	h5py, numpy	<pre>real_data = f[path][0][0] ppm = (frequencies / spectrometer_freq)[-1]</pre>
2. Data Processing	Noise reduction and smoothing	scipy.signal, scipy.ndimage	<pre>denoised_data = medfilt(real_data, kernel_size=5) smoothed_data = gaussian_filter(real_data, sigma=2)</pre>
3. Range Analysis	Calculate integral values within a specified ppm range	numpy	<pre>mask = (ppm >= r[0]) & (ppm <= r[1]) sum(data_to_analyze[mask])</pre>
4. Graph Generation	Visualize data	matplotlib.pyplot	<pre>plt.plot(ppm, real_data) plt.savefig("original_plot.png")</pre>
5. Noise Evaluation	Compare noise levels using standard deviation	numpy	<pre>np.std(original_data) noise_reduction = (...)</pre>
6. Statistical Analysis	Compute mean, max, and min values	numpy	<pre>np.mean(data) np.max(data)</pre>
7. PDF Report Generation	Compile results and graphs into a report	reportlab.platypus	<pre>doc = SimpleDocTemplate(...) flowables.append(Image(...))</pre>

Example: Python-Based Data Analysis Using JASON and HDF5

This section presents a simple workflow for NMR data analysis and report generation using Python, leveraging the structure of HDF5 files:

1. Load NMR data
2. Process and analyze the data
3. Automatically generate graphs and reports

Details of each step and the Python libraries used are summarized in Table 1.

This example serves as a practical reference for understanding HDF5 file structures and performing spectral analysis.

Reading NMR Data and Calculating the ppm Axis

The function shown in Figure 1 demonstrates how to extract real signal data and compute the chemical shift (ppm axis) from NMR measurement files in .j5[4] format.

This approach is applicable to various spectra, including ¹³C NMR, and can be adapted as a general-purpose data loading function by accounting for differences in file structure and metadata attributes.

Example of Path Specification

HDF5 files are structured like hierarchical databases, consisting of folders and datasets organized in layers.

Below is an example of how to specify a path:

```
path = "JasonDocument/NMR/NMRData/0/DataPoints"
```

```
real_data = f[path]['0'][(0)] # Real component of the signal
```

Note:

The string "JasonDocument/NMR/NMRData/0/DataPoints" represents the hierarchical structure within the file.

The expression f[path]['0'][(0)] accesses the data associated with the key '0' under the specified path.

Metadata such as SpecInfo can also be accessed using similar path specifications.

Example of Generated Report

Figure 2 shows a sample report generated from the workflow described in Table 1.

For more details on libraries such as NumPy, SciPy, and ReportLab, please refer to their respective user guides.

```
import h5py
import numpy as np

def extract_nmr_data(jjh5_file):
    with h5py.File(jjh5_file, 'r') as f:
        # Retrieve real part of the spectral data
        # → Extract the "real" component of the NMR signal
        path = "JasonDocument/NMR/NMRData/0/DataPoints"
        real_data = f[path]['0'][(0)] # Real part

        # Retrieve measurement conditions (metadata) to calculate the ppm axis
        # → Necessary information to construct the horizontal axis (ppm) of
        # the spectrum
        specinfo_path = "JasonDocument/NMR/NMRData/0/SpecInfo"
        attrs = f[specinfo_path].attrs

        # Sweep Width (SW): Frequency range measured (Hz)
        sweep_width = attrs["SW"][0]

        # Spectrometer Frequency (MHz): Central frequency of the instrument
        spectrometer_freq = attrs["SpectrometerFrequencies"][0]

        # Spectrum Reference Position (Hz)
        spectrum_ref = attrs["SpectrumRef"][0]

        # Calculate frequency axis
        num_points = len(real_data)
        start_freq = spectrum_ref - (sweep_width / 2)
        stop_freq = spectrum_ref + (sweep_width / 2)
        frequencies = np.linspace(start_freq, stop_freq, num_points)

        # Convert to ppm axis (divide by MHz), then reverse to increase from
        # left to right
        ppm = (frequencies / spectrometer_freq)[::-1]

    return ppm, real_data
```

Figure 1. Excerpt from the NMR Data Extraction Script



Figure 2. Example of Report Output

2. Use Case of the External Command Function

Overview of the External Command Function

The External Command function in JASON allows users to execute external Python scripts directly from the Processing tab in the JASON GUI.

This feature is particularly useful for automation and custom data processing.

Here, we introduce an example of using Python scripts from the downloadable package External NMR Processing Scripts, available on the official JASON website.

Sample Code Access

Sample scripts can be downloaded from the following JASON official resource page:

<https://www.jeoljason.com/resources-external-nmr-processing-scripts/>

Figure 3 shows a list of Python scripts included in the External NMR Processing Scripts directory.




Figure 3. External NMR Processing Scripts

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

```

from argparse import ArgumentParser

# Prepare to handle command-line arguments
# + Allows users to specify filenames or parameters when running the script
parser = ArgumentParser(description=
    "A tool for editing parameters in JASON data files.")

# Target filename (JSON file in HDF5 format)
parser.add_argument("-f", "--filename", action="store", default="$TMPFILE",
    help="Path to the target JASON HDF5 file")

# Name of the parameter to edit (e.g., SpectrometerFrequencies)
parser.add_argument("-p", "--par", action="store",
    default="SpectrometerFrequencies",
    help="Name of the parameter to edit")

# Dimension index to edit (0-7)
# + For multidimensional data, specifies which dimension to modify
parser.add_argument("-d", "--dim", action="store", type=int, default=0,
    help="Index of the dimension to edit (0-7)")

# New value (float)
# + The new value to assign to the specified parameter
parser.add_argument("-v", "--vNum", action="store", type=float, default=500.0,
    help="New value to set for the parameter")

# Parse the arguments and store them in a variable
args = parser.parse_args()

```

Figure 4. Excerpt from the `jsonParEdit.py` Script

About the `jsonParEdit.py` Script

`jsonParEdit.py` is a Python script designed to edit parameters within JASON data files, such as SW, SpectrumRef, and SpectrometerFrequencies.

It can be executed from the GUI to modify values interactively.

Figure 4 shows a portion of the `jsonParEdit.py` source code.

This script allows users to freely modify the parameters mentioned above.

Usage instructions are included at the top of the script and should be reviewed before execution.

How to Use

The operation procedure is summarized in Table 2.

An example of the GUI settings and input format is shown in Figure 5.

Table 2. Operation Steps for Using External Command in JASON GUI

Step	Process Description
1	Add the External Command option to the Processing tab in the JASON GUI.
2	In the Cmd: field of the External Command tab, specify the full path to the Python executable.
3	In the Arguments: field, specify the full path to the Python script to be executed (in this case, <code>jsonParEdit.py</code>).
4	Define the processing parameters to be executed by <code>jsonParEdit.py</code> . For example, to change the sweep width (SW) to 5000 Hz, use the following command: C:\Users\externalNMRProcessing\python\jsonParEdit.py -f \$TMPFILE -p SW -d 0 -v 5000.0
5	After confirming the settings in the Processing tab, click Apply to execute the command.




Figure 5. Example Configuration for `jsonParEdit.py`




Figure 6. Example of SW Output from
jasonParEdit.py




Figure 7. Example of SpectrumRef Output
from jasonParEdit.py

Output Example

Figure 6 presents an example of the spectral output after processing.

The updated spectrum is reflected in the JASON GUI and can be compared with the original data.

Error Messages

If incorrect settings are provided during script execution, clear error messages are displayed to help identify the issue:

- Unsupported parameter name: Error: parameter: XXX is not supported by this script.
- Nonexistent attribute: Error: parameter: XXX does not exist.
- Invalid dimension index: Error: dimension index must be in range of 0-7.

These messages ensure that any issues encountered after clicking the Apply button can be quickly diagnosed, making the script safe and user-friendly.

Example: Adjusting SpectrumRef with jasonParEdit.py

Figure 7 illustrates an example of adjusting the SpectrumRef parameter using jasonParEdit.py.

The operation method follows the same steps as described in Table 2.

Below is an example of a full path argument used to set the observation center to 0.0 ppm:

```
C:\Users\externalNMRProcessing\python\jasonParEdit.py -f
$TMPFILE -p SpectrumRef -d 0 -v -0.0
```

Supplement: Exploring HDF5 File Structure

If the structure of an HDF5 file is unclear, the following types of exploration scripts can be used to inspect the types and locations of stored data:

Structure Exploration Script Example: Allows listing of all datasets within the file.

```
import h5py

with h5py.File("your_file.jjh5", "r") as f:
    f.visit(print)
```

Metadata Inspection: Enables retrieval of measurement conditions and other metadata (e.g., frequency, sweep width).

```
with h5py.File("your_file.jjh5", "r") as f:
    for key in f["JasonDocument/NMR/NMRData/0/SpecInfo"].attrs:
        print(key, f["JasonDocument/NMR/NMRData/0/SpecInfo"].attrs[key])
```

Tips

For interactive exploration and inspection of HDF5 file contents, users can use HDFView, a free graphical browser provided by The HDF Group.

HDFView allows users to navigate the hierarchical structure of HDF5 files, view datasets, and examine metadata without any programming.

Summary

By combining JASON, HDF5, and Python, efficient processing and analysis of NMR data becomes possible.

- The hierarchical structure of HDF5 enables flexible data extraction.
- The External Command function facilitates seamless integration with external programs such as Python.
- Custom processing and report generation enhance analytical flexibility.
- Understanding HDF5 structure and utilizing exploration scripts are key to advanced applications.

References

JASON Official Website: <https://www.jeoljason.com/>

External NMR Processing Scripts: <https://www.jeoljason.com/resources-external-nmr-processing-scripts/>

Python Documentation: <https://docs.python.org/>

h5py Library: <https://www.h5py.org/>

NumPy Library: <https://numpy.org/>

ReportLab Library: <https://www.reportlab.com/documentation/>

[1] JEOL Analytical Software Network

[2] Python is a trademark or registered trademark of the Python Software Foundation.

[3] SpecInfo is a structure that stores metadata related to NMR measurements, such as experimental conditions. SW refers to the sweep width (measurement range), SpectrumRef indicates the center frequency, and SpectrometerFrequencies represents the observation frequency.

[4] .jjh5 is a file format used in JASON, based on the HDF5 structure.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.