


Airtight measurement tool for diffusion NMR of electrolytes

Product used : Nuclear Magnetic Resonance (NMR)

Introducing the NMR measurement tool for anaerobic samples such as battery electrolytes. Samples can be easily packed even in a glove box, where sampling operations are restricted. It is compatible with 5mm-diameter solution NMR sample cups, making it easy to perform NMR measurements in an airtight conditions. Recommended for diffusion measurements of solution electrolytes and solid electrolytes.



- ※ A torque screwdriver is required separately.
- ※ Please contact our sales for information on compatible NMR sample tubes.